Product Recommendation
of COVID-19 Neutralizing Antibody

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
Corona Virus (COVID-19) Neutralizing
Antibody Detection Kit (Colloidal Gold)

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The immune system will produce corresponding antibodies after being stimulated by the virus. Some antibodies can bind to the neutralizing epitope of the virus particles, making the virus lose the ability to bind to the receptor, preventing the virus from infecting cells, so as to "neutralize" the virus. These antibodies are called neutralizing antibodies.

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The detection of neutralizing antibody of novel coronavirus can reflect the individual's recovery and the effect of vaccine, but the traditional virus neutralization experiment needs to use live virus and cells, and has high requirements for laboratory safety and operator skills. This time-consuming and demanding traditional neutralization detection method does not apply to COVID-19 with fast transmission speed and wide coverage. Therefore, the kit provides a simple, time-consuming and convenient auxiliary diagnostic method for the detection of neutralizing antibody against neocoronavirus.

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

1. 检测实验必须在 18~30°C 之间进行。
The test must be done at temperature between 18 ~ 30°C.
2. 检测结果只能表明该样本中的中和抗体的水平，不能准确量化中和抗体滴度。
The test results can only show the level of neutralizing antibody in the sample, and can not accurately quantify the neutralizing antibody titer.
3. 快速测试应密封并保存在干燥的地方，并应在打开包装后立即使用。
The rapid test should be sealed and kept in dry place and should be used as soon as the packing is opened.
4. 快速检测结果仅供临床参考，不应成为临床诊断和治疗的唯一依据。
The results of rapid test are only for clinical reference and should not be the only basis for clinical diagnosis and treatment.
5. 废弃样本和检测后的试剂卡应作为潜在的传染源处理。
Waste samples and test should be treated as potential infectious agents.
6. 质控线出现时间不作为判断检测线结果的依据，应在 15–20 分钟内判读结果。
The result should be read strictly within a time limit of 15–20 minutes.
7. 该试剂只能用于检测接种疫苗后或感染恢复后的中和抗体水平，不能用于评估抗病毒保护的有效性。
This reagent can only be used to detect the level of neutralizing antibody after vaccination or recovery from infection, and cannot be used to evaluate the effectiveness of antiviral protection.

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Beijing Zhongjian Antai Diagnostic Tech. Co., Ltd.

Dirección(*)

No.29 Qingfeng Road West, Daxing biological medicine base

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Nombre Comercial (*)

CORONA VIRUS (COVID-19) NEUTRALIZING ANTIBODY DETECTION KIT (COLLOIDAL GOLD)

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CE EC Declaration of Conformity CE

Manufacturer: Beijing Zhongjian Antai Diagnostic Technology Co., Ltd
29 Qingfeng Road West,Daxing biological medicine base,Zhongguancun Science and Technology Park, Daxing, Beijing 102629,P.R.China

Whose Single Authorized EU- Representative: Riomavix S.L.
Add.: Calle de Almansa 55, 1D, Madrid 28039 Spain

Product Name: Corona Virus (COVID-19) Neutralizing Antibody Detection Kit (Colloidal Gold)

Classification : **Others of ANNEX II of IVDD**

Conformity Assessment Route: **Annex III**

We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.

General applicable directives:

In Vitro Diagnostic Medical Devices DIRECTIVE 98/79/EC

Harmonized standards:

EN ISO 13485:2016,EN ISO 15223-1:2016,EN ISO 14971:2012, EN 13641: 2002,
EN ISO 18113-1:2011, EN ISO 18113-2:2011, EN 13612: 2002, EN ISO 23640:2015

Signature:

ZHAO ZHENGCHUN

Mr. Zhao Zhengchun
Deputy General manager of Operations
China, July 6,2021

Name:

Mr. Zhao Zhengchun

Title:

Deputy General manager of Operations

Place/Date:

China, July 6,2021

Declaration of Conformity

Manufacturer: Beijing Zhongjian Antai Diagnostic Technology Co., Ltd
29 Qingfeng Road West,Daxing biological medicine base,Zhongguancun Science and Technology Park, Daxing, Beijing 102629,P.R.China

Whose Single Authorized EU- Representative: Riomavix S.L.
Add.: Calle de Almansa 55, 1D, Madrid 28039 Spain

Product Name: Corona Virus (COVID-19) Neutralizing Antibody Detection Kit (Colloidal Gold)

Classification : **Others of ANNEX II of IVDD**

Conformity Assessment Route: **Annex III**

We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.

General applicable directives:

In Vitro Diagnostic Medical Devices DIRECTIVE 98/79/EC

Harmonized standards:

EN ISO 13485:2016,EN ISO 15223-1:2016,EN ISO 14971:2012, EN 13641: 2002,
EN ISO 18113-1:2011, EN ISO 18113-2:2011, EN 13612: 2002, EN ISO 23640:2015

Signature:

Name:

Mr. Zhao Zhengchun

Title:

Deputy General manager of Operations

Place/Date:

China, July 6,2021

RIO MAVIX

NOTIFICATION OF REGISTRATION

This is to certify that, according to the European Council Directive 98/79/EC, Riomavix S.L. performed all notification duties and responsibilities as the European Authorized Representative:

MANUFACTURER: Beijing Zhongjian Antai Diagnostic Technology Co., Ltd

ADDRESS: No.29 Qingfeng Road West, Daxing biological medicine base, Zhongguancun Science and Technology Park, Daxing, Beijing 102629, P.R.China

The manufacturer has provided Riomavix S.L. with all the appropriate declaration according to the European Council Directive 98/79/EC including the Declaration of Conformity confirming that its in vitro diagnostic medical device, as stipulated here below, is fulfilling the essential requirements of the European Council Directive 98/79/EC.

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Classification: Others

Where the manufacturer affix the CE mark to the device listed they must ensure that all the essential requirements of European Council Directive 98/79/EC are met.

The notification of abovementioned device has been completed by the European Authorized Representative in Spain. The Spain Competent Authority is notified of the manufacture's device and has allocated registration. The registration number is RPS/1632/2021

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Beijing Zhongjian Antai Diagnostic Technology Co., Ltd.

No.29 Qingfeng Road, Daxing Bio-pharmaceutical Industry Base, Zhongguancun Science & Technology Park, Daxing District, Beijing, 102629, P.R.China

Index of Symbol

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Riomavix S.L.

Add.: Calle de Almansa 55, 1D, Madrid 28039 Spain

E-mail: leis@riomavix.com

Tel.: +34 658 396 230

Version Number : 1.1

Effective Date: June 4, 2021

Corona Virus (COVID-19) Neutralizing Antibody Detection Kit (Colloidal Gold)

【Product name】

Corona Virus (COVID-19) Neutralizing Antibody Detection Kit (Colloidal Gold)

【Packing specification】

20 tests / box

【Intended use】

This product is for professional use and is used for in vitro detection of antibody levels in human blood samples after being vaccinated with the Novel Coronavirus vaccine or recovered from the Novel Coronavirus. But it cannot be used to evaluate the effectiveness of anti-virus protection.

【Test principle】

This kit adopts the principle of colloidal gold immunochromatography technology. The detection area (T line) of the nitrocellulose membrane is coated with mouse anti-human IgG monoclonal antibody, and the quality control area (C line) is coated with Mouse anti-RBD monoclonal antibody. Colloidal gold-labeled COVID-19 recombinant antigen is coated on a gold-labeled pad. When testing the sample, the COVID-19 neutralizing antibody in the sample combines with the COVID-19 Ag labeled with colloidal gold (Au) to form (Au- COVID-19Ag-[COVID-19-Neutralizing antibody]) immune complex. As the chromatographic complex flows forward inside the nitrocellulose membrane, it combines with the coated mouse anti-human IgG monoclonal antibody when passing through the detection area to form "(Au- COVID-19Ag-[COVID-19- Neutralizing antibody])- [mouse anti-human IgG monoclonal antibody]" form agglutination, and the remaining colloidal gold labeled COVID-19 recombinant antigen combines with the mouse anti-RBD monoclonal antibody coated at the quality control line to agglutinate Color rendering. For example, if the sample does not contain COVID-19 neutralizing antibodies, which prevents the formation of immune complexes in the detection area, only the quality control area will show color.

【Principal component】

The test consists of a membrane strip coated with mouse anti-human IgG on the test line, and the dye pad contains colloidal gold conjugated with COVID-19 recombinant antigen.

1. Test cassette (20 tests)
2. Sample diluent (4mL × 1bottle)

【Storage conditions and expiry date】

- 4~30°C dry, keep away from light, valid for 12 months.
- The product should be stored in dry condition under 4~30°C and kept away from light. When humidity is less than 60%, it should be used within 1 hour and when humidity is more than 60%, it should be used immediately.
- Expiration date and lot number are shown in the label.

【Sample requirements】

- The whole blood should be venous blood. Serum samples can be collected by vein in a conventional manner.

The plasma samples can be treated with heparin, sodium citrate and EDTA. Above samples can be placed under 2~8°C for 5 days. Samples under -20°C can be stored for at least 3 months.

- Samples should avoid hemolysis or repeated freezing-thawing. If the sample is turbid or has precipitation, it should be centrifuged or filtered to clarify before testing.

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This package insert must be read completely before performing the test. Please restore the reagent and sample to room temperature before inspection. Experimental humidity should be less than 60%, the experiment temperature is 18~30°C. Test procedure is as follows:

1. Remove the test card from the aluminum foil bag, mark the sample and put it on the horizontal work table.
2. Take 20 μL serum, plasma samples or 40 μL whole blood to be directly added to the add hole or the bottom of the indicator arrow, and then add the sample diluent 100 μL (about 2-3 drops).
3. Result should be read at 15-20 minutes, negative results must be confirmed at the end of 20 minutes.

【Interpetation of assay result】

1. Positive: color development of the detection line and quality control line. It is suggested that the neutralizing antibody of the novel coronavirus (COVID-19) detected is positive.
2. Negative: Only the quality control line develops color in the detection window. It is suggested that the concentration of neutralizing antibody for novel coronavirus (COVID-19) in the sample has not reached the level of detection.
3. Invalid results: The control line has no red stripe. The invalid results should resume experiment again, and the test should be in strict accordance with the instructions operation, if the test result is still invalid, please contact local suppliers or customer service with our company for technical consultation.

【Limitations of the test method】

1. The product is only used for the detection of whole blood, serum or plasma.
2. This product inspection result is only for clinical reference, should not serve as the only basis for clinical diagnosis and treatment.
3. This reagent can only be used to detect the neutralizing antibody titer after vaccination or recovery from infection, and cannot be used to evaluate the effectiveness of antiviral protection.
4. The results of sample testing are related to factors such as sample collection, testing, transportation and storage. Any error will affect the accuracy of the results.
5. Interfering substances: 1) When the bilirubin concentration ≤ 20mg/dL, the hemoglobin content ≤ 500mg/dL, and the triglyceride content ≤ 1500mg/dL, it will not interfere with the test results of this product; 2) when the antinuclear antibody titer ≤ 1:320, rheumatism factor ≤ 500IU/mL will not interfere with the test results of this product.

【Product performance indicator】

1. Compliance rate of positive quality control products: the positive internal quality control compliance rate should be 8/8.
2. Compliance rate of negative quality control products: the negative internal quality control compliance rate should be 8/8.
3. Minimum detection limit: the minimum detection limit of quality control products should not be lower than 1:8. (NT50 ≥ 10)
4. Repeatability: test 2 internal repeatability quality control products, each test for 10 times, results should be

positive.

【Cautions】

1. The operation shouldn't be done if the condition of kit and samples are not restored to 18 ~ 30°C, in case it will affect the accuracy of the results.
2. The positive samples obtained by the rapid test should be confirmed by other methods.
3. The rapid test should be sealed and kept in dry place. The test bar should be tested as soon as possible after being removed from the packaging, so as to avoid placing it in the air for too long, causing the damp.
4. The deepness of the test line color is not necessarily associated with the titer of the antibody in the sample, and the results of the interpretation after 20 minutes are invalid.
5. The results of rapid test are only for clinical reference and should not be the only basis for clinical diagnosis and treatment.
6. Waste samples and test should be treated as potential infectious agents.
7. The appearing time of the Control Line should not be taken as the time basis for judging the results of test line. The color rendering results should be observed and judged within a time limit of 15-20 minutes.
8. The rapid test is only used for in vitro diagnosis.

Beijing Zhongjian Antai Diagnostic Technology Co., Ltd.

No.29 Qingfeng Road, Daxing Bio-pharmaceutical Industry Base, Zhongguancun Science & Technology Park, Daxing District, Beijing, 102629, P.R.China

Authorized representative in the European Community

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
The management system of

Beijing Zhongjian Antai
Diagnostic Technology Co., Ltd.

No. 29 Qingfeng West Road, Daxing Bio-pharmaceutical Industry Base,
Zhongguancun Science & Technology Park, Daxing District,
Beijing City, 102609, P.R. China

has been assessed and certified as meeting the requirements of

ISO 13485:2016
EN ISO 13485:2016

For the following activities

Design, Development and Manufacture of In Vitro Diagnostic Kits based on Gold immunochromatography Assay and Fluorescence Immunochromatography Assay; Design, Development and Manufacture of Fluorescence Immunochromatography Analyzer

This certificate is valid from 18 January 2021 until 17 January 2024
and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 29 November 2023
Issue 3. Certified since 18 January 2018

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
Beijing Zhongjian Antai Diagnostic Technology Co., Ltd. is a controlling member company of Beijing Hinspire Healthcare Group Co., Ltd. Hinspire Healthcare Group are involved in high-tech fields such as biomedicine, medical instruments, diagnostic reagents, mobile medical treatment, intellectual endowment, high polymer materials and so on.

Hinspire Healthcare Group relies on Internet plus medical innovation incubation and Beijing international science and technology innovation base, create a healthy industry cluster; Establish investment and financing platform through the incubation industry capital link.

Beijing Zhongjian Antai Diagnostic Technology Co., Ltd. is based on the research and development, registration, production and sales of in-vitro diagnostic reagents. At present, there are two technical platforms of colloidal gold and fluorescent immunity, the colloidal gold platform is mainly type III products, while the fluorescent immune platform is mainly type II quantitative products. The company has 12 utility model patents, 1 invention patent, 1 design patent and 12 software copyrights.. As a high-tech enterprise with GMP standard production plant, the company belongs to the national high-tech enterprises and high-tech enterprises in Zhongguancun, had a full set of license.

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SGS United Kingdom Ltd