

DOA Cup Reader

Provide an intuitive user interface and high performance to scan the results on Round Test Cups.

Key features:

- Stand alone device with its own touch screen and keyboard
- Compact size and light weight
- Quick boot — no lamp warm up time
- Cup scan < 40 seconds
- Easy operation with minimized steps
- Auto-recognition of drug test cup type and number of strips
- Provides “positive”, “negative” and “invalid” judgment to each test strip
- Highly sensitive to weak lines, more objective than eye readout
- Data saved as .csv and .pdf format
- Export data through USB thumb drive or PC-link
- Brief on-station record review
- Future upgrade options to more types of rapid diagnostic cups and flat cards



EASY



FAST



ACCURATE



- 2nd floor, No. 3, 8th Street, Qiantang Area, 310018 Hangzhou
- Web: www.prometheusbio.com
- Tel: 86-571-86908151
- Fax: 86-571-86908139
- E-mail: sales@prometheusbio.com

Other configurations may also be available. Please call us for more information.




TEST NAME	CALLIBRATOR	CUT-OFF
6-MAM	6-Monoacetylmorphine	10ng/ml
K3	AB-PINACA	300ng/ml
UR144	Oxazepam	300ng/ml
CAF	Buprenorphine	10ng/ml
GBPT	Benzoyllecgonine	300ng/ml
MDPV	Delta-9-THC-COOH	50ng/ml
Zolpidem	D-Metamphetamine	1000/500ng/ml
CLO	Methadone	300ng/ml
METH	MDMA	500ng/ml
Meprobamate	Morphine	300ng/ml
HDCD	Morphine	2000ng/ml

TEST NAME	CALLIBRATOR	CUTOFF	CE
Amphetamine	D-Amphetamine	1000ng/ml	*
Barbiturates	Socobarbital	300ng/ml	*
Benzodiazepine	Oxazepam	300ng/ml	*
Buprenorphine	Buprenorphine	10ng/ml	*
Cocaine	Benzoyllecgonine	300ng/ml	*
Marijuana	Delta-9-THC-COOH	50ng/ml	
Metamphetamine	D-Metamphetamine	1000/500ng/ml	
Methadone	Methadone	300ng/ml	*
MDMA	MDMA	500ng/ml	*
Morphine	Morphine	300ng/ml	*
Opiates 2000	Morphine	2000ng/ml	*
Oxycodone	Oxycodone	100ng/ml	*
Phencyclidine	Phencyclidine	25ng/ml	*
Propoxyphene	Propoxyphene	300ng/ml	*
Tricyclic Antidepressants	nortriptyline	1000ng/ml	*
Alcohol/ALC	Alcohol	0.04%B.A.C	
Fentanyl/FEN	Fentanyl	100ng/ml	
K2/Spice	JWH-018/JWH-073	50/25 ng/ml	*
EDDP	EDDP	300ng/ml	*
Tramadol	Tramadol	300ng/ml	*
Methaqualone	Methaqualone	300ng/ml	
MDPV	MDPV	300ng/ml	
Zolpidem	Zolpidem	10ng/ml	
Ketamine	Ketamine	1000ng/ml	
Cotinine	Cotinine	200ng/ml	

Name
John Smith

Donor ID
123456

1	Negative
2	Positive
3	Negative
4	Negative
5	Negative
6	Positive
7	Negative
8	Negative
9	Negative
10	Negative
11	Negative
12	Negative



Other configurations may also be available. Please call us for more information.