

# DRUG OF ABUSE

OUR PRODUCTS

One Step Rapid Test for Drug of Abuse  
Easy Fast Accurate



- 2nd floor, No. 3, 8th Street, Qiantang Area, 310018 Hangzhou, Zhejiang, China
- Office: 86-571-86908151
- Fax: 86-571-86908139
- Email: [sales@prometheusbio.com](mailto:sales@prometheusbio.com)
- Website: [www.prometheusbio.com](http://www.prometheusbio.com)

## Product Introduction:


### Excellent Design:

- Most cost-effective drug test cup with built-in adulteration strip
- Simple procedure with no additional steps required to activate test
- Screw lid for non-leakage
- Various product format including strip, cassette, dipcard and cup
- Both single and multiple drug tests are available
- Both urine and saliva tests are available
- Customized configurations available

### Reliable Results:

- Easy-to-read result display in strong color codes
- Long-lasting results stable up to 1h
- 99% accuracy proved by clinical studies
- 24 months shelf life
- CE approved
- Set to SAMHSA cut-off levels

**Please find the common configurations as following.  
For more information please feel free to contact with us.**



Negative (-) : both line with C and T



Positive (+) : only show C line



Invalid: without C line

TEST NAME	CALIBRATOR	CUTOFF
* Amphetamine	D-Amphetamine	1000ng/mL
* Ketamine	Ketamine	1000ng/mL
* Barbiturates	Secobarbital	300ng/mL
* Benzodiazepines	Oxazepam	300ng/mL
* Buprenorphine	Buprenorphine	10ng/mL
* Cocaine	Benzoyllecgonine	300ng/mL
* Marijuana	Delta-9-THC-COOH	50ng/mL
* Methamphetamine	D-Methamphetamine	1000ng/mL
* Methadone	Methadone	300ng/mL
* MDMA	MDMA	500ng/mL
* Morphine	Morphine	300ng/mL
* Opiates 2000	Morphine	2000ng/mL
* Oxycodone	Oxycodone	100ng/mL
* Phencyclidine	Phencyclidine	25ng/mL
* Propoxyphene	Propoxyphene	300ng/mL
* Tricyclic Antidepressants	Nortriptyline	1000ng/mL
* Cotinine	Cotinine	200ng/mL
* Fentanyl/FEN	Fentanyl	100ng/mL
* K2/Spice	JWH-018/JWH-073	50ng/mL 25ng/mL
* EDDP	EDDP	300ng/mL
* Tramadol	Tramadol	300ng/mL
APAP	Acetaminophen	2000ng/mL
...	...	...

\* CE marked

### 1、Urine Test



### 2、Saliva Test



Other configurations may also be available. Please call us for more information.