



Drug Screening Test Cup with Adulteration Test

Following urine collection, peel the label and view the test result in minutes.

Excellent Design

- Patent-pending drug screening test cup with built-in adulteration strip
- Simple procedure with no additional steps required to activate test
- Screw lid for non-leakage
- Test format from 1 to 17 drugs per cup
- Customized configurations available

Reliable Results

- Easy-to-read result displayed in strong color codes
- Long-lasting results stable up to 1 hour
- 99% accuracy proved by clinical studies
- 24 months shelf life
- Set to SAMHSA cut-off levels
- FDA approved
- **CLIA Waived**

TEST NAME	CALIBRATOR	CUTOFF
Δ Amphetamine	D-Amphetamine	1000 ng/mL
Δ Barbiturates	Secobarbital	300 ng/mL
Δ Benzodiazepines	Oxazepam	300 ng/mL
Δ Buprenorphine	Buprenorphine	10 ng/mL
Δ Cocaine	Benzoylcegonine	300 ng/mL
Δ Marijuana	Delta-9-THC-COOH	50 ng/mL
Δ Methadone	Methadone	300 ng/mL
Δ Methamphetamine	D-Methamphetamine	1000 ng/mL
Δ MDMA	MDMA	500 ng/mL
Δ Opiates 300	Morphine	300 ng/mL
Δ Opiates 2000	Morphine	2000 ng/mL
Δ Oxycodone	Oxycodone	100 ng/mL
Δ Phencyclidine	Phencyclidine	25 ng/mL
Δ Propoxyphene	Propoxyphene	300 mg/mL
Δ Tricyclic Antidepressants	Nortriptyline	1000 ng/mL
Δ Amphetamine	D-Amphetamine	300 ng/mL

TEST NAME	CALIBRATOR	CUTOFF
Δ Cocaine	Benzoylcegonine	150 ng/mL
Δ Methamphetamine	D-Methamphetamine	500 ng/mL
+ Alcohol/ALC	Alcohol	0.04%
+ APAP	Acetaminophen	2000 ng/mL
+ Cotinine/COT	Cotinine	200 ng/mL
+ EDDP	EDDP	300 ng/mL
+ Ethyl Glucuronide	Ethyl Glucuronide	500 ng/mL
+ Fentanyl/FEN	Fentanyl	100 ng/mL
+ K2/Spice	JWH-018/JWH-073	50 ng/mL
+ Ketamine	Ketamine	1000 ng/mL
+ Lysergic Acid Diethylamide	Lysergic Acid Diethylamide	10 ng/mL
+ MDPV	MDPV	300 ng/mL
+ Methaqualone	Methaqualone	300 ng/mL
+ Ritalin	Ritalin	300 ng/mL
+ Tramadol	Tramadol	300 ng/mL
+ Zolpidem	Zolpidem	10 ng/mL

Δ 510K

* CLIA waived

+ For forensic use only

Other configurations may also be available. Please call us for more information