



Drug Test Cup Reader



The UCP-SR100 scanner provides an intuitive user interface and high performance to scan the results on UCP's Drug Test Cups (Round Cups)

Feature

- Stand alone device with its own touch screen and keyboard
- Compact size and light weight
- Quick boot – no lamp warm up time
- Cup scan < 40 seconds
- Easy operation with minimized steps
- Auto-recognition of drug test cup type and number of strips
- Provides “positive”, “negative” and “invalid” judgement to each test strip
- Highly sensitive to weak lines, more objective than eye readout
- Data saved as .csv and .pdf format
- Export data through USB thumb drive or PC-link
- Brief on-station record review
- Future upgrade options to more types of rapid diagnostic cups and flat cards.

2nd floor, No. 3, 8th Street, Qiantang Area, 310018 Hangzhou, Zhejiang, China
T: +86(0571)86908151 sales@prometheusbio.com

Prometheus™
Your Trusted Rapid Tests