

# Infectious Disease

OUR PRODUCTS

One Step Rapid Test  
for Infectious Disease

Easy Fast Accurate

**Prometheus**<sup>TM</sup>  
Your Trusted Rapid Tests

Web: [www.prometheusbio.com](http://www.prometheusbio.com)


E-mail: [sales@prometheusbio.com](mailto:sales@prometheusbio.com)

Address: 2nd floor, No. 3, 8th Street, Qiantang Area,  
310018 Hangzhou, Zhejiang, China

Other configurations may also be available. Please call us for more information.

## Product Introduction:

- Cost-effective method with less turn around time.
- CE marked products
- Easy-to-read result display in strong color codes
- For various kind of disease



Positive (+) : both line with C and T



Negative (-): only show C line



Invalid: without C line

## What We Provide?



DESCRIPTION	ITEMS	SPECIMEN	TEST/KIT	REGISTRATION
Respiratory Pathogens	Influenza A/B Antigen	Secretion	1, 10, 25, 40	*
	Respiratory syncytial Virus Antigen	Secretion	1, 10, 25, 40	*
	Group A Streptococcus Antigen	Secretion	1, 10, 25, 40	*
	M.tuberculosis Antibody	WB/S/P	1, 10, 25, 40	*
	M.tuberculosis IgG/IgM	WB/S/P	1, 10, 25, 40	*
	M.pneumoniae IgM	WB/S/P	1, 10, 25, 40	*
	C.pneumoniae IgM	WB/S/P	1, 10, 25, 40	*
	Respiratory virus Antigen Combined-Test 3 in 1	Secretion	1, 20, 40	*
	Respiratory virus Antigen Combined-Test 4 in 1	Secretion	1, 20, 40	*
Vector Borne Pathogens	Malaria P.f/P.v	WB	1, 10, 25, 40	*
	Malaria P.f/Pan	WB	1, 10, 25, 40	*
	Dengue IgG/IgM	WB/S/P	1, 10, 25, 40	*
	Dengue NS1	WB/S/P	1, 10, 25, 40	*
	Dengue NS1+IgG/IgM	WB/S/P	1, 10, 25, 40	*
Digestive System Pathogens	S.typhi Antigen	Feces	1, 10, 25, 40	*
	H.pylori Antibody	WB/S/P	1, 10, 25, 40	*
	H.pylori Antigen	Feces	1, 10, 25, 40	*
	Hepatitis A virus IgM	WB/S/P	1, 10, 25, 40	*
	Adenovirus Antigen	Feces	1, 10, 25, 40	*
Maternal and child health	Rotavirus Antigen	Feces	1, 10, 25, 40	*
	C.trachomatis IgM	WB/S/P	1, 10, 25, 40	*
	C.trachomatis Ag	Secretion	1, 10, 25, 40	*
	Herpes simplex Virus IgG/IgM	WB/S/P	1, 10, 25, 40	*
Other	Group B Streptococcus Antigen	Secretion	1, 10, 25, 40	*
	Human Brucella IgM/IgG	WB/S/P	1, 10, 25, 40	*
	C. difficile toxicin A/B Antigen	Feces	1, 10, 25, 40	*
	Fecal Occult Blood	Feces	1, 25, 40 10, 25, 50, 100	*
	Monkeypox Virus Antigen	Secretion	1, 10, 25, 40	*

\* CE Marked

Other configurations also be available. Please call us for more information.