



PrometheusTM
「Your Trusted Rapid Tests」

Add: 2nd floor, No. 3, 8th Street, Qiantang Area, 310018 Hangzhou, Zhejiang, China

Web: www.prometheusbio.com

Tel: 86- 571- 86908151

Fax: 86- 57186908139

E-mail: sales@prometheusbio.com

PrometheusTM
「Your Trusted Rapid Tests」

PRODUCT CATALOG





PROMETHEUS

COMPANY PROFILE

With decades of committed service to the rapid test, Prometheus is committed to a true partnership with customers by providing easy, fast and accurate diagnostics products. Focusing on technology, innovation and quality control, our manufacturing covers a board range of areas: Drug Test Series; Women's Health Series; Cardiac Markers Series; Infectious Disease Series and Tumor Markers Series. The DOA Cup Reader, scanning data from the urine cup, can lead the drug testing market, and the Up-Converting Fluorescence Quantitative Meter-spot measurement also has a considerable advantage in the market. Guided by our precise technology, we are continuously exchanging knowledge with our customers, partners, and clinicians to address the critical issues of improving care and lowering cost.

Prometheus plant was granted ISO 9001 :2008 & ISO 13485:2003 certification. In November 2012, Prometheus' production facility was inspected by US FDA and no 483 form was issued. High standards and compliance, to which Prometheus faithfully adheres, allow our customers to have the utmost confidence in the validity of their project.



CONTENTS

- 01 DRUG OF ABUSE
- 07 INFECTIOUS DISEASES
- 09 WOMEN'S HEALTH
- 10 CARDIAC MARKERS
- 11 TUMOR MARKERS
- 13 NEW PRODUCT

DRUG OF ABUSE /Urine Test




Product Introduction

Excellent Design:

- Most cost-effective drug test cup with built-in adulteration strip
- Simple procedure required to activate test
- Screw lid for non-leakage
- Various product format including strip, cassette, dipcard and cup
- Urine, saliva and hair tests are available
- Customized configurations available

Reliable Results:

- Easy-to-read result display in strong color codes
- Over 99% accuracy proved by clinical studies
- 24months shelf life
- CE, FDA and NMPA approved
- Set to SAMHSA cut-off levels



Description	Specimen	Format	Sensitivity	Test /Kit	Registration	
Amphetamine/AMP	Urine	Strip	1000 ng/mL	50	△	*
		Cassette	300 ng/mL	25	△	*
			1000 ng/mL	25	△	*
		Card	300 ng/mL	25	△	*
			1000 ng/mL	25	△	*
			300 ng/mL	25	△	*
N -acetyl-p-aminophenol/APAP	Urine	Strip	2000 ng/mL	50		
		Cassette	2000 ng/mL	25		
		Card	2000 ng/mL	25		
Barbiturates/BAR	Urine	Strip	300 ng/mL	50	△	*
		Cassette	300 ng/mL	25	△	*
		Card	300 ng/mL	25	△	*
Benzodiazepines/BZO	Urine	Strip	300 ng/mL	50	△	*
		Cassette	300 ng/mL	25	△	*
		Card	300 ng/mL	25	△	*
Buprenorphine/BUP	Urine	Strip	10 ng/mL	50	△	*
		Cassette	10 ng/mL	25	△	*
		Card	10 ng/mL	25	△	*
Cocaine/COC	Urine	Strip	300 ng/mL	50	△	*
			150 ng/mL		△	
			100 ng/mL			*
		Cassette	300 ng/mL	25	△	
			150 ng/mL			
			100 ng/mL			*
		Card	300 ng/mL	25	△	*
			150 ng/mL		△	
			100 ng/mL			
Cotinine/COT	Urine	Strip	200 ng/mL	50		*
		Cassette	200 ng/mL	25		*
		Card	200 ng/mL	25		*
EDDP	Urine	Strip	300 ng/mL	50		*
			100 ng/mL			*
			300 ng/mL	25		*
		Cassette	100 ng/mL	25		*
			300 ng/mL			*
			100 ng/mL			*
Fentanyl /FEN	Urine	Strip	100 ng/mL	50		
		Cassette	100 ng/mL	25		*
		Card	100 ng/mL	25		
Ketamine/KET	Urine	Strip	1000 ng/mL	50	✓	*
		Cassette	1000 ng/mL	25	✓	*
		Card	1000 ng/mL	25	✓	*

△: FDA Approved *: CE Approved ✓: CFDA Approved

Description	Specimen	Format	Sensitivity	Test /Kit	Registration		
K2	Urine	Strip	50 ng/mL	50			*
			25 ng/mL				*
		Cassette	50 ng/mL	25			*
			25 ng/mL				*
		Card	50 ng/mL	25			*
			25 ng/mL				*
LSD	Urine	Strip	10 ng/mL	50			
			3 ng/mL				
		Cassette	10 ng/mL	25			
			3 ng/mL				
		Card	10 ng/mL	25			
			3 ng/mL				
Methadone/MTD	Urine	Strip	300 ng/mL	50	△		*
			200 ng/mL				
			50 ng/mL				
		Cassette	300 ng/mL	25	△		*
			200 ng/mL				
			50 ng/mL				
		Card	300 ng/mL	25	△		*
			200 ng/mL				
			50 ng/mL				
Methamphetamines/MET/mAMP	Urine	Strip	1000 ng/mL	50	△	✓	*
			500 ng/mL		△		
			300 ng/mL				
		Cassette	1000 ng/mL	25	△	✓	*
			500 ng/mL		△		
			300 ng/mL				
		Card	1000 ng/mL	25	△	✓	*
			500 ng/mL		△		
			300 ng/mL				
Methylenedioymethamphetamine/ MDMA	Urine	Strip	500 ng/mL	50	△	✓	*
		Cassette	500 ng/mL	25	△	✓	*
		Card	500 ng/mL	25	△	✓	*
Morphine/MOP	Urine	Strip	300 ng/mL	50	△	✓	*
		Cassette	300 ng/mL	25	△	✓	*
		Card	300 ng/mL	25	△	✓	*
Methaqualone/MQL	Urine	Strip	300 ng/mL	50			
		Cassette	300 ng/mL	25			
		Card	300 ng/mL	25			
Methylphenidate/MTHP	Urine	Strip	300 ng/mL	50			
		Cassette	300 ng/mL	25			
		Card	300 ng/mL	25			

Description	Specimen	Format	Sensitivity	Test /Kit	Registration							
Opiates/OPI	Urine	Strip	2000 ng/mL	50	△		*					
			200 ng/mL									
			100 ng/mL									
			50 ng/mL									
		Cassette	2000 ng/mL	25	△		*					
			200 ng/mL									
			100 ng/mL									
			50 ng/mL									
		Card	2000 ng/mL	25	△		*					
			200 ng/mL									
			100 ng/mL									
			50 ng/mL									
Oxycodone/OXY	Urine	Strip	300 ng/mL	50	△		*					
			100 ng/mL									
		Cassette	300 ng/mL	25	△		*					
			100 ng/mL									
		Card	300 ng/mL	25	△		*					
			100 ng/mL									
		Marijuana/THC	Urine	Strip	50 ng/mL	50	△	✓	*			
					40 ng/mL							
					20 ng/mL							
				Cassette	50 ng/mL	25	△	✓	*			
					40 ng/mL							
					20 ng/mL							
Card	50 ng/mL			25	△	✓	*					
	40 ng/mL											
	20 ng/mL											
Phencyclidine/PCP	Urine			Strip	25 ng/mL	50	△		*			
				Cassette	25 ng/mL					25	△	*
				Card	25 ng/mL							
Propoxyphene/PPX	Urine	Strip	300 ng/mL	50	△		*					
		Cassette	300 ng/mL					25	△	*		
		Card	300 ng/mL									
Tricyclic Antidepressants/TCA	Urine	Strip	1000 ng/mL	50	△							
		Cassette	1000 ng/mL					25	△	*		
		Card	1000 ng/mL									
Tramadol/TRA	Urine	Strip	300 ng/mL	50			*					
		Cassette	300 ng/mL					25		*		
		Card	300 ng/mL									
Alcohol/ALC	Urine	Strip	0.04%	50								
		Cassette	0.04%					25				
		Card	0.04%									
Ethyl Glucuronide/EtG	Urine	Strip	500 ng/mL	25								

Description	Specimen	Format	Sensitivity	Test /Kit	Registration	
Multi-Drug Screen Card	Urine	2-14 Drugs	See Insert	25	△	*
		2- 12 Drugs with S.V.T.	See Insert	25		*
MultiDrug Screen Cassette	Urine	2-12 Drugs	See Insert	25	△	*
MultiDrug Screen Cup	Urine	2-12 Drugs	See Insert	25	△	*
		2-10 Drugs with S.V.T.	See Insert	25		
		2- 15 Drugs	See Insert	25	△	*
		2-15 Drugs with S.V.T.	See Insert	25		
		2- 18 Drugs	See Insert	25	△	*
		2-16 Drugs with S.V.T.	See Insert	25		

△: FDA Approved *: CE Approved

▶

ADULTERATION

Description	Specimen	Format	Test /Kit
pH+S.G.+Ox	Urine	Strip	50
NIT+Glu+CRE	Urine	Strip	50
pH+S.G.+Ox+CRE	Urine	Strip	50
S.G.+CRE	Urine	Strip	50
S.G+CRE+pH	Urine	Strip	50
CRE+pH	Urine	Strip	50
Creatinine/CRE	Urine	Strip	50
NIT+S.G+CRE+pH	Urine	Strip	50
NIT+OX+CRE+pH	Urine	Strip	50
NIT+OX+S.G.+pH	Urine	Strip	50
OX+S.G.+CRE	Urine	Strip	50

DRUG OF ABUSE /Saliva Test



Description	Specimen	Format	Sensitivity	Test /Kit	Registration
Alcohol /ALC	Saliva	Strip	> 0.02% B.A.C	50	*
		Cassette	> 0.02% B.A.C	25	*
Amphetamine/ AMP	Saliva	Cassette	50 ng/mL	25	*
Barbiturates/BAR	Saliva	Cassette	50 ng/mL	25	*
Buprenorphine/BUP	Saliva	Cassette	10 ng/mL	25	*
Cocaine/COC	Saliva	Cassette	20 ng/mL	25	*
Methamphetamines/MET/ mAMP	Saliva	Cassette	50 ng/mL	25	*
Methadone/MTD	Saliva	Cassette	35 ng/mL	25	*
Methylenedioxymethamphetamine/ MDMA	Saliva	Cassette	60 ng/mL	25	*
Opiates /OPI	Saliva	Cassette	40 ng/mL	25	*
Oxycodone/OXY	Saliva	Cassette	50 ng/mL	25	*
Propoxyphene/PPX	Saliva	Cassette	50 ng/mL	25	*
Phencyclidine/PCP	Saliva	Cassette	10 ng/mL	25	*
Marijuana/THC	Saliva	Cassette	12 ng/mL	25	*
Ketamine/KET	Saliva	Cassette	100 ng/mL	25	*
Cotinine/COT	Saliva	Cassette	40 ng/mL	25	*
MultiDrug Screen Cassette	Saliva	2-12Drugs	See Insert	25	*
MultiDrug Screen Cup	Saliva	2-12Drugs	See Insert	25	*

*: CE Approved


INFECTIOUS DISEASES



▶ Product Introduction

Excellent Design:

- Cost- effective method with less turn around time
- Easy-to-read result display in strong color codes
- For various kind of diseases
- CE marked products



Description	Specimen	Format	Catalog	Test /Kit	Registration
Chlamydia	Swab/Urine	Cassette	ICH-502	25	
Gonorrhea	Swab	Cassette	ICH-502	25	*
Strep A	Throat Swab	Strip	ST-501	25	*
		Cassette	IST-502	25	*
HBsAg	S/P	Strip	IHBsg-301	50	
		Cassette	IHBsg-302	25	
	WB/S/P	Strip	IHBsg-401	50	
		Cassette	IHBsg- -402	25	
HBsAb	S/P	Strip	IHBsb- 301	50	
		Cassette	IHBsb--302	25	
HIV 1/2/O	WB/S/P	Strip	IHI401	50	
		Cassette	IHI--402	25	
HIV 1/2	WB/S/P	Strip	IHI-A401	50	
		Cassette	IHI-A402	25	
HCV Antibody	S/P	Strip	IHC -301	50	
		Cassette	IHC -302	25	
HCV Antibody	WB/S/P	Strip	IHI- A401	50	
		Cassette	IHI- A402	25	
Treponema Pallidum	WB/S/P	Strip	ITP-401	50	*
		Cassette	ITP-402	25	*
	S/P	Cassette	ITP -302	25	*
Tuberculosis Mononucleosis	WB/S/P	Cassette	ITB-402	25	
	WB/S/P	Cassette	IMO- 402	25	
Malaria P.f/Pan	WB	Strip	IMA-401	50	*
		Cassette	IMA-402	25	*
Malarla P.f	WB	Strip	IMA-A401	50	
		Cassette	IMA -A402	25	
H.pylori Antigen	Feces	Cassette	IHP-602	25	
H.pylori Antibody	S/P	Cassette	IHP-302	25	*
H.pylori Antibody	WB/S/P	Strip	IHP 401	50	*
		Cassotte	IHP 402	25	*
Influenza A& B	Swab	Strlp	IFL-501	25	*
		Cassotte	IFL-502	25	*
Rotavirus	Feces	Cassette	IAD -602	25	*
Adenovirus	Feces	Cassette	IRO-602	25	*
Rotavlrus & Adenovirus	Feces	Cassette	IRA-625	25	
Dengue IgG/IgM	WB/S/P	Cassette	IDE- 402	25	*
2019-nCoV IgG/IgM	WB/S/P	Cassette	ICOV-402	25	
2019-nCoV IgG/IgM Test Cassette	WB/S/P	Cassette	ICOV-402	25	* NEW
2019-nCoV Antigen Test Cassette	Swab	Cassette	ICOV-502	25	* NEW
2019-nCoV Neutraliza- tion Antibody Test Cassette	WB/S/P	Cassette	INCOV-402	25	* NEW
2019-nCoV&Flu A/B Combo Antigen Test Cassette	Swab	Cassette	ICOVF-502	25	* NEW

*: CE Approved


For Distribution Outside US Only

WOMEN'S HEALTH



Product Introduction

- 25 mIU is available
- Designed for OTC testing
- Reliable result in 5 minutes, simple to read
- CE marked
- Easy-to-read result display in strong color codes
- Upcoming product: digital pregnancy test with weeks indicator



Description	Specimen	Format	Sensitivity	Registration
HCG Pregnancy	Urine	Strip	25 mIU/mL	△ ✓ *
		Cassette	10 mIU/mL	△ ✓ *
		Card	25 mIU/mL	△ ✓ *
		Midstream	10 mIU/mL	△ ✓ *
	Urine/Serum	Strip	25 mIU/mL	△ ✓ *
		Cassette	10 mIU/mL	△ ✓ *
		Card	25 mIU/mL	△ ✓ *
		Midstream	10 mIU/mL	△ ✓ *
LH Ovulation	Urine	Strip	25 mIU/mL	✓
		Cassette	25 mIU/mL	✓
		Card	25 mIU/mL	✓
		Midstream	25 mIU/mL	✓


△: FDA Approved *: CE Approved ✓: CFDA Approved For Distribution Outside US Only

CARDIAC MARKERS



Product Introduction

- Easy-to-read result display in strong color codes
- Reliable result in 10 - 20 minutes, simple to read
- Long-lasting results stable up to 0.5h
- CE marked products
- 18 months shelf life



Description	Specimen	Format	Catalog	Test/Kit	Registration
Cardiac Troponin I (cTnI)	WB/S/P	Cassette	CTI-402	25	*
Myoglobin/Creatine Kinase MB/ Cardiac Troponin I Combo	WB/S/P	Cassette	CMA-435	25	*

*: CE Approved


For Distribution Outside US Only

TUMOR MARKERS



Product Introduction

- Easy-to-read result display in strong color codes
- Reliable result in 5 minutes, simple to read
- Long-lasting results stable up to 1h
- CE marked products
- 24 months shelf life



Description	Specimen	Format	Catalog	Sensitivity	Test/Kit	Registration
Fecal Occult Blood FOB	Feces	Strip	TFO-601	50 ng/mL	25	*
				6 Ug/g	25	
	Feces	Cassette	TFO-602	50 ng/mL	25	*
				6 Ug/g	25	

*: CE Approved

For Distribution Outside US Only

UPCOMING NEW PRODUCT

Description	Specimen	Sensitivity	Callibrator
6-MAM	Urine	10 ng/mL	6-Monoacetylmorphine
K3	Urine	300 ng/mL	AB-PINACA
K4	Urine	300 ng/mL	UR-144- 5-pentanoic acid
CAF	Urine	8000 ng/mL	Caffeine
GBPT	Urine	1000 ng/mL	Gabapentin
MDPV	Urine	300 ng/mL	Methylenedioxypropyvalerone
ZOLP	Urine	10 ng/mL	Zolpidem
CLO	Urine	300 ng/mL	Clonazepam
METH	Urine	500 ng/mL	S(-)-Methcathinone Hcl
Meprobamate	Urine	1000 ng/mL	Meprobamate
HDCD	Urine	100 ng/mL	Hydrocodone
CLZP	Urine	300 ng/mL	Klonopin
PGB	Urine	2000 ng/mL	Pregabalin

NEW PRODUCT
DOA CUP READER



EASY



FAST



ACCURATE

▶ KEY FEATURES

- Stand alone device with its own touch screen and keyboard
- Compact size and light weight
- Quick boot - no lamp warm up time
- Cup scan < 40 seconds
- Easy operation with minimized steps
- Auto-recognition of drug test cup type and number of strips
- Provides " positive", "negative" and "invalid" judgment to each test strip
- Highly sensitive to weak lines, more objective than eye readout
- Data saved as .CSV and .pdf format
- Export data through USB thumb drive or PC-link
- Brief on- -station record review
- Future upgrade options to more types of rapid diagnostic cups and flat cards



PROMETHEUS SOLO



▶ KEY FEATURES

- Novel fluorescence nanoparticles have virtually invincible illumination and zero crosstalk to sample matrix in signaling
- Patented particle chemistry, surface chemistry and formulation chemistry leading to silky lateral flow with very little particle residuals on nitrocellulosemembrane
- High sensitivity, high precision, fully quantitative detection over broad range of biomarker concentration
- Ultra-fast detection, all tests are done within 3 -12 minutes for seamless diagnostic decision making process
- Test compatible with various bodily fluids such as serum, plasma and wholeblood

YOUR TRUSTED RAPID TESTS